



Recombinant tRNA-specific adenosine deaminase (tadA)

Product Code	CSB-MP733385SMZ
Abbreviation	tadA
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5XE14
Product Type	Recombinant Protein
Immunogen Species	Streptococcus pyogenes serotype M6 (strain ATCC BAA-946 / MGAS10394)
Purity	≥85% (SDS-PAGE)
Sequence	MPYSLEEQTY FMQEALKEAE KSLQKAEIPI GCVIVKDGEI IGRGHNAREE SNQAIMHAEM MAINEANAHE GNWRLDITTL FVTIEPCVMC SGAIGLARIP HVIYGASNQK FGGADSLYQI LTDERLNHRV QVERGLLAAD CANIMQTFRR QGRERKKIAK HLIKEQSDPF D
Source	Mammalian cell
Target Names	tadA
Protein Names	Recommended name: tRNA-specific adenosine deaminase EC= 3.5.4.-
Expression Region	1-171
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.