



Recombinant Rat Transgelin-2 (TagIn2)

Product Code	CSB-BP733395RA
Abbreviation	TagIn2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5XFX0
Product Type	Recombinant Protein
Immunogen Species	Rattus norvegicus (Rat)
Purity	≥85% (SDS-PAGE)
Sequence	ANRGPSTYGL SREVQQKIEK QYDPDLEQIL IQWITTQCRK GVSQPQPGRE NFQNWLDKGT VLCELINSLY PEGQAPVKKI QASTMAFKQM EQISQFLQAA ERYGINTTDI FQTVDLWEGK NMACVQRTLML NLGGLAVARD DGLFSGDPNW FPKKSSENPR NFSDNQLQEG KNVIGLQMGT NRGASQAGMT GYGMPRQIL
Source	Baculovirus
Target Names	TagIn2
Protein Names	Recommended name: Transgelin-2
Expression Region	2-199
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Target Details	This protein is a homolog of the protein transgelin, which is one of the earliest markers of differentiated smooth muscle. The function of this protein has not yet been determined.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.