



# Recombinant *Oryza sativa* subsp. japonica Delta-aminolevulinic acid dehydratase, chloroplastic (HEMB)

<b>Product Code</b>	CSB-EP733510OFG
<b>Abbreviation</b>	HEMB
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q5Z8V9
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	<i>Oryza sativa</i> subsp. japonica (Rice)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	VSSEQ EAAPAVRAPS GRTIEECEAD AVAGRFPAPP PLVRPKAPEG TPQIRPLDLT KRPRRNRSP ALRAAFQETT ISPANLVLPL FIHEGEDDAP IGAMPGCYRL GWRHGLLDEV YKSRDVGVNS FVLFPKVPDA LKSQSGDEAY NDNGLVPRTI RLLKDKFPDI VVYTDVALDP YSSDGHGIV REDGVIMNDE TVYQLCKQAV SQARAGADV VSPDMMDGRV GAIRAALDAE GFHDVSIMSY TAKYASSFYG PFREALDSNP RFGDKKTYQM NPANYREALL ETAADAEAGA DILLVKPGLP YLDVIRLLRD NSALPIAAYQ VSGEYSMIKA GGALNMIDEE KVMMESLMCL RRAGADIILT YFARQAANVL CGMRSN
<b>Source</b>	E.coli
<b>Target Names</b>	HEMB
<b>Protein Names</b>	Recommended name: Delta-aminolevulinic acid dehydratase, chloroplastic Short name= ALADH EC= 4.2.1.24 Alternative name(s): Porphobilinogen synthase
<b>Expression Region</b>	46-426
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.