



# Recombinant Chicken BRCA1-A complex subunit Abraxas (FAM175A)

<b>Product Code</b>	CSB-YP733517CH
<b>Abbreviation</b>	FAM175A
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q5ZHS0
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Gallus gallus (Chicken)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MEGESTSALL SGFVFGALAF QHLSTDSLTE GFLLDGDKGE AKNSITDSQM DDVEVVYTID IQKHIPCYQL FSFYNAGEL NEPALKKILS GRKKSIVIGWY KFRNRDQIM TFRERLLHKN LQSHLSNQGL VFLLLTSSVM TESCSTYRLE HALHRPQEGL FQKVPLVVTN LGMAEQQGYR TVSGSCASSG FVRAVKQHRS EFFYEDGSLQ EVHKINEMYA TLQEELKMC SDVEVSERSV EKLLEVSQLE KEEINRKKQH KISSGGNKDQ LEEPKENVLL CQALRTFFPS SDLQTCIVSF KGRISKNCCK IDHNINIMDK LTLMVEERDF TEAETKLVTK RKVRVTTTGP KSLKKLRSLE LDQELHQDEE DCNQETKLAL SSAETDEEAL ENPKDTNEYS YSPTF
<b>Source</b>	Yeast
<b>Target Names</b>	ABRAXAS1
<b>Protein Names</b>	Recommended name: BRCA1-A complex subunit Abraxas Alternative name(s): Coiled-coil domain-containing protein 98 Protein FAM175A
<b>Expression Region</b>	1-405
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.