



Recombinant Chicken Programmed cell death protein 2-like (PDCD2L)

Product Code	CSB-MP733524CH
Abbreviation	PDCD2L
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5ZID2
Product Type	Recombinant Protein
Immunogen Species	Gallus gallus (Chicken)
Purity	≥85% (SDS-PAGE)
Sequence	MAAAAPPVLR GLRDAAMVGP CRGGGEPPGW ATNKLGGGRAD ALPSVRPGCP RCGACGAALA HLVQVYCP LG ASPFHRLANV FACAESACWG RPQSWKVLRS QSLEARGQDA LGCGSKQKEE CNFSAKDWCA EADDWGGCDG TESSTCALQM LGLNAGMSSS ASEEACVSVQ LQELHLSEAA DVSGSLNTDL LLREEMVMAT SAPVFHSHYYI SVVDEADYAG FLDTHANKL LKEYQQREGV DLEHLMSESF AGEDSNEKYE KTKVRS GDHT FHKFMKRISV CPEQILRYSW GGQPLFITCP PANFDNSIPA CSNCGSNRVF EFQLMPALVS MLWGDADLSV EFGTVIVYTC ERSCWPTNQQ TPLEEFIFVQ EDPDQKLFK
Source	Mammalian cell
Target Names	PDCD2L
Protein Names	Recommended name: Programmed cell death protein 2-like
Expression Region	1-379
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.