



Recombinant Chicken Protein FAM172A (FAM172A)

Product Code	CSB-EP733544CH
Abbreviation	FAM172A
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5ZK44
Product Type	Recombinant Protein
Immunogen Species	Gallus gallus (Chicken)
Purity	≥85% (SDS-PAGE)
Sequence	VPDEKEK TTALKDLLSR IDLDELMKKD EPPLEFPDTL EGFEYSFNEK GQLRHAKTGE PFVFNYREDL HRWNQKRYEA LGEIITKHVY GLLKECHLK KVTLPVDATE NEPKSFIFMS EDALTNPKL LVLHNGV RAGQWARRLI INEDLDSGTQ IPYIKKATEE GYGVIVLNP ENYIEVEKTK AQVQLSSDSS DEPAEKREK DKIQKETKKR RDFYEKYRNP QKEKETMQMY IRDNGSPEEH AIYVWDHFIS QSAAENVFFV AHSYGGLAFV ELMIQRETEV KNKVTAVALT DSVHNVWHQE AGKTIREWMR ENCCNWVSSS EPLDTSVESM LPDCPRVSAG TERHELTSWK SFPSIFKFFS EAMEAKNSSV KPAPTRRSNR IKYEEL
Source	E.coli
Target Names	FAM172A
Protein Names	Recommended name: Protein FAM172A
Expression Region	24-416
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.