



Recombinant Chicken Protein Muted homolog (MUTED)

Product Code	CSB-BP733545CH
Abbreviation	MUTED
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5ZK77
Product Type	Recombinant Protein
Immunogen Species	Gallus gallus (Chicken)
Purity	≥85% (SDS-PAGE)
Sequence	MSGAGPASPG RSGSALPLPL PLSLSLTGDR KREVPAAGTL LQPIIKDVGE IYSRLLDHRP VIQGEIRYFV KEFEEKRGLR ELRVLENLKK TIFETNEDIL PKCKQSMHDN LNEVLQRLQA ANDTIHKLQE KEHEKRKIQA DKLMAGEEKR IAQQETFLQE QQSKREEVDE EHRKAMERLK EQYSEMEKEL AKYVSF
Source	Baculovirus
Target Names	BLOC1S5
Protein Names	Recommended name: Protein Muted homolog
Expression Region	1-196
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Target Details	This gene encodes a component of BLOC-1 (biogenesis of lysosome-related organelles complex 1). Components of this complex are involved in the biogenesis of organelles such as melanosomes and platelet-dense granules. A mouse model for Hermansky-Pudlak Syndrome is mutated in the murine version of this gene. Some transcripts of the downstream gene TXNDC5 overlap this gene, but they do not contain an open reading frame for this gene.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.