



Recombinant Chicken Charged multivesicular body protein 6 (CHMP6)

Product Code	CSB-YP733558CH
Abbreviation	CHMP6
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5ZL55
Product Type	Recombinant Protein
Immunogen Species	Gallus gallus (Chicken)
Purity	>85% (SDS-PAGE)
Sequence	GNLFGKRKR SRVTEQDKAV LQLKQQRDKL RQYQKRISLG LERERELARQ LLKEGKKEKA MLLLKKKRYQ EQLLDKTDNQ ISNLERMVQD IEFTQIEMKV IEGLKIGNEC LNKMHQVMSI EEVERIIGET QDAVEYQRQI DEILAGSLTE EDEDAILLEEL NAITQEQLLEL PDVPSEPLPE EPPEATPVKN RPKPELVAAS
Source	Yeast
Target Names	CHMP6
Protein Names	Recommended name: Charged multivesicular body protein 6 Alternative name(s): Chromatin-modifying protein 6
Expression Region	2-200
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.