



Recombinant *Xanthomonas campestris* pv. *campestris* Protein RecA (recA)

Product Code	CSB-YP733636XAY
Abbreviation	recA
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q60101
Product Type	Recombinant Protein
Immunogen Species	<i>Xanthomonas campestris</i> pv. <i>campestris</i> (strain ATCC 33913 / DSM 3586 / NCPPB 528 / LMG 568 / P 25)
Purity	≥85% (SDS-PAGE)
Sequence	MDENKKRALS AALSQIEKQF GKGSVMRMGD RVIEAVEVIP TGSLMLDIAL GIGGLPKGRV VEIYGPESSG KTTTLQAIA ECQKLGGTAA FIDAEHALDP IYAAKLG VNV DDLLSQPDT GEQALEIADM LVRSSSVDIV VIDSVAAALTP KAEIEGEMGD QLPGLQARLM SQALRKLTGN IKRSNTLVVF INQLRMKIGV MMPGQSPEVT TGGNALKFYA SVRLDIRRIG AIKKGDEIIG NQTKIKVVKN KLAPPFKQVI TEILYGEGIS REGELIDMGV EAKLVDKAGA WYSYGDERIG QGKDNARGYL RDN PQVAIKL EAE LREKFQP AEAPREAGET ESE
Source	Yeast
Target Names	recA
Protein Names	Recommended name: Protein RecA Alternative name(s): Recombinase A
Expression Region	1-343
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.