



Recombinant *Methylococcus capsulatus* Enolase 2 (eno2)

Product Code	CSB-EP733675MFM
Abbreviation	eno2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q604M4
Product Type	Recombinant Protein
Immunogen Species	<i>Methylococcus capsulatus</i> (strain ATCC 33009 / NCIMB 11132 / Bath)
Purity	≥85% (SDS-PAGE)
Sequence	MSKIVDILAR EILDSRGNPT VQAEVILDSG AEGSAMVPSG ASTGAREAIE LRDGDASRYG GKGVLKAVEN VRGPIKAALT GMDAADQAAI DRRLIELDGS DNKGVLGANA ILAVSLATAH AAAADAKKPL YAYLNRSGEF LLPVPMMNII NGGAHADNSI DMQEFMILPV GAPSREALR YGAEVFHALK KVLSDRGLAT GVGDEGGFAP NLPSEAAIG IILEAIEKAG YRPGEDICLG LDVASSEFYS DGIYTLASEG RRFTSEEFSD HLAAWVDKYP IVSIEDGMAE NDWHGWGIHT DRLGRRIQLV GDDLFTNPA ILKQGIEARI ANSILIKVNQ IGTLTETLEA IHIARKAGYT SVVSHRSGET EGTTIADIAV ATCTGQIKTG SLSRSDRIAK YNRLKIEEE LGDKARYGGR GAVKNLA
Source	E.coli
Target Names	eno2
Protein Names	Recommended name: Enolase 2 EC= 4.2.1.11 Alternative name(s): 2-phospho-D-glycerate hydro-lyase 2 2-phosphoglycerate dehydratase 2
Expression Region	1-427
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.