



# Recombinant Methylococcus capsulatus Enolase 2 (eno2)

|                          |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                          |
|--------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Code</b>      | CSB-YP733675MFM                                                                                                                                                                                                                                                                                                                                                                                                                                                                                          |
| <b>Abbreviation</b>      | eno2                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                     |
| <b>Storage</b>           | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.                                                                                                                                                                                                                       |
| <b>Uniprot No.</b>       | Q604M4                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                   |
| <b>Product Type</b>      | Recombinant Protein                                                                                                                                                                                                                                                                                                                                                                                                                                                                                      |
| <b>Immunogen Species</b> | Methylococcus capsulatus (strain ATCC 33009 / NCIMB 11132 / Bath)                                                                                                                                                                                                                                                                                                                                                                                                                                        |
| <b>Purity</b>            | ≥85% (SDS-PAGE)                                                                                                                                                                                                                                                                                                                                                                                                                                                                                          |
| <b>Sequence</b>          | MSKIVDILAR EILDSRGNPT VQAEVILDSG AEGSAMVPSG ASTGAREAIE<br>LRDGDASRYG GKGVLKAVEN VRGPIKAALT GMDAADQAAI DRRLIELDGS<br>DNKGVLGANA ILAVSLATAH AAAADAKKPL YAYLNRSGEF LLPVPMMNII<br>NGGAHADNSI DMQEFMILPV GAPSREALR YGAEVFHALK KVLSDRGLAT<br>GVGDEGGFAP NLPSEAAIG ILEAIEKAG YRPGEDICLG LDVASSEFYS<br>DGIYTLASEG RRFTSEEFSD HLAAWVDKYP IVSIEDGMAE NDWHGWGIHT<br>DRLGRRIQLV GDDLFTNPA ILKQGIEARI ANSILIKVNQ IGTLTETLEA<br>IHIARKAGYT SVVSHRSGET EGTTIADIAV ATCTGQIKTG SLRSRDRIAK<br>YNRLKIEEE LGDKARYGGR GAVKNLA |
| <b>Source</b>            | Yeast                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                    |
| <b>Target Names</b>      | eno2                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                     |
| <b>Protein Names</b>     | Recommended name: Enolase 2 EC= 4.2.1.11 Alternative name(s): 2-phospho-D-glycerate hydro-lyase 2 2-phosphoglycerate dehydratase 2                                                                                                                                                                                                                                                                                                                                                                       |
| <b>Expression Region</b> | 1-427                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                    |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.                                                                                                                                                                                                                                                                                                                                                                                                      |
| <b>Tag Info</b>          | Tag type will be determined during the manufacturing process.                                                                                                                                                                                                                                                                                                                                                                                                                                            |
| <b>Protein Length</b>    | full length protein                                                                                                                                                                                                                                                                                                                                                                                                                                                                                      |
| <b>Reconstitution</b>    | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.                                                                                              |
| <b>Shelf Life</b>        | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life                                                                                                                                                                                                                                                                        |



of lyophilized form is 12 months at -20°C/-80°C.