



Recombinant Methylococcus capsulatus

Tryptophan synthase alpha chain (trpA)

Product Code	CSB-MP733676MFM
Abbreviation	trpA
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q604P4
Product Type	Recombinant Protein
Immunogen Species	Methylococcus capsulatus (strain ATCC 33009 / NCIMB 11132 / Bath)
Purity	≥85% (SDS-PAGE)
Sequence	MSRLTSTFQV LRHANRKALI PFITAGDPEP GFTVPVMHAM VEAGADVIEL GVPFSDPMAD GPVIQKASER ALVHHMSLGK TLDLVAEFRR ADRTTPIVLM GYLNPVECMG YQRFVDRAKD AGVDGVLTV D MPPEEAGEFV PMLHRAGLDP IFLLAPNSSD ARVQRMAALG RGYLYYVSLK GVTGASHLDL AEVESRIAH I RSLTDLP IGV GFGVKNAETA AAIGRFADAV VVGSALVEKI AAAQDRPEQA RKDIVDLLAV MRRALDA
Source	Mammalian cell
Target Names	trpA
Protein Names	Recommended name: Tryptophan synthase alpha chain EC= 4.2.1.20
Expression Region	1-267
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.