



Recombinant Methylococcus capsulatus RNA 3'-terminal phosphate cyclase (rtcA)

Product Code	CSB-MP733729MFM
Abbreviation	rtcA
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q609M1
Product Type	Recombinant Protein
Immunogen Species	Methylococcus capsulatus (strain ATCC 33009 / NCIMB 11132 / Bath)
Purity	≥85% (SDS-PAGE)
Sequence	MTLLRIDGSY GEGGGQIIRT ALSLAALTGT AVRLENIRAR RRPGLAAQH LTAVKAVALL CDARIGGAEL GSQTLEFIGP RPVRAGDYTL DVGEAREGGS AGAVMLVLQT VLPPLALAEAG ASTVSLRGGT HVPMSPPFDY VREVWLPTLA RMGVRAELSL VRSGWYPVGK GEVRLVVAGG RRSLEALKAR HRGTLQAVTG RALAANLPAH IPERMAARAR AVLADAGIAA AIEPAVLQAA CAGAGLFLTA RYDHCLAGFT ALGRRGKSAE RVAEEACGAL LAHHRSGAAL DQHLADQIVL PAALCREESL FSVERITPHL TTNAWVDRF GLARVDVRLA ECGTGLVRIH GNPAGVCHAH P
Source	Mammalian cell
Target Names	rtcA
Protein Names	Recommended name: RNA 3'-terminal phosphate cyclase Short name= RNA cyclase Short name= RNA-3'-phosphate cyclase EC= 6.5.1.4
Expression Region	1-361
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.