



Recombinant *Macaca fascicularis* Torsin-1A (TOR1A)

Product Code	CSB-EP733754MOV-B
Abbreviation	TOR1A
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q60HG2
Product Type	Recombinant Protein
Immunogen Species	<i>Macaca fascicularis</i> (Crab-eating macaque) (<i>Cynomolgus</i> monkey)
Purity	≥85% (SDS-PAGE)
Sequence	VEPISLGLAL AGVLTGYIYP RLYCLFAECC GQKRSL SREA LQKDLDNKLF GQHLAKKIIL NAVFGFINNP KPKKPLTSL HGWTGTGKNF VSKIIAENIY EGGLNSDYVH LRVATLHFPN ASNITLYKDQ LQLWIRGNVS ACARSIFIFD EMDKMHAGLI DAIKPF LDYY DLVDGVS YQK AIFIFLSNAG AERITDVALD FWRSGKQRED IKLKDIEHAL SVSVFNKNS GFWHSSLIDR NLIDYFVPFL PLEYKHLKMC IRVEMQSRGY ETNEDIVSRV AEEMTFFPKE ERVFSKDGCK TVFTKLDYYY DD
Source	<i>E.coli</i>
Target Names	TOR1A
Protein Names	Recommended name: Torsin-1A Alternative name(s): Dystonia 1 protein Torsin family 1 member A
Expression Region	21-332
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.