



Recombinant Putative H/ACA ribonucleoprotein complex subunit 2-like protein (CBG18231)

Product Code	CSB-BP733763CXX
Abbreviation	CBG18231
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q60YI3
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis briggsae
Purity	≥85% (SDS-PAGE)
Sequence	MGKRNLDETM NESAVSEIAE DVAALTTEKD EYQALCELVN PIAQPLANRK LAKKVYKLIK KAAAGEKTLR EGIKDVQKEL RKNEKGICIL AGNVSPIDVY SHIPAICEEK EIPYVYIPSR EQLGLAVGHR RPSILIFVKP SADFAQELYDE VAETLHHLTV DAA
Source	Baculovirus
Target Names	CBG18231
Protein Names	Recommended name: Putative H/ACA ribonucleoprotein complex subunit 2-like protein
Expression Region	1-163
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.