



Recombinant Mouse Neurexophilin-2 (Nxph2)

Product Code	CSB-EP733777MO
Abbreviation	Nxph2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q61199
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	≥85% (SDS-PAGE)
Sequence	EEVIHNTSVDWEDRTVP ETLVGNLFHS RITSPLRLFV KQPPDPKPSY ADNTKNFWDW LANITEIQEQ LARTKRRPIV KTGKFKKMFG WGDHFSNIKT VKLNLLITGK IVDHGNGTFS VYFRHNSTGL GNVSVSLVPP SKVVEFEISP QSTLETKESK SFNCHIEYEK TDRAKKTALC NFDPSKICYQ EQTQSHVSWL CSKPFKVICI HIIFYSVDYK LVQKVCPDYN YHSETPYLSF G
Source	E.coli
Target Names	Nxph2
Protein Names	Recommended name: Neurexophilin-2
Expression Region	23-261
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.