



Recombinant Rat Egl nine homolog 3 (Egln3)

Product Code	CSB-BP733831RA
Abbreviation	Egln3
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q62630
Product Type	Recombinant Protein
Immunogen Species	Rattus norvegicus (Rat)
Purity	>85% (SDS-PAGE)
Sequence	MPLGHIMRLD LEKIALEYIV PCLHEVGFCY LDNFLGEVVG DCVLERVKQL HYNGALRDGQ LAGPRAGVSK RHLRGDQITW IGGNEEGCEA INFLLSLIDR LVLYCGSRLG KYYVKERSKA MVACYPGNGT GYVRHVDNPN GDGRCITCIY YLKNWDAKL HGGVLRIFPE GKSFVADVEP IFDRLLFSWS DRRNPHEVQP SYATRYAMTV WYFDAEERAE AKKKFRNLTR KTESALAKD
Source	Baculovirus
Target Names	Egln3
Protein Names	Recommended name: Egl nine homolog 3 EC= 1.14.11.29 Alternative name(s): Hypoxia-inducible factor prolyl hydroxylase 3 Short name= HIF-PH3 Short name= HIF-prolyl hydroxylase 3 Short name= HPH-3 Prolyl hydroxylase do
Expression Region	1-239
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.