



Recombinant RNA-binding protein pno-1 (pno-1)

Product Code	CSB-EP733836CXX
Abbreviation	pno-1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q626C1
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis briggsae
Purity	≥85% (SDS-PAGE)
Sequence	MATSSSAFDD EFPMEEGMPE LLDDDEDVPST LPSLLEQNMD TALKSNEFKL VKRKRRAGNA VDVVMEDISQ PDDSTTDSPD ADAEQKPTKR SKGSKGESRV VPVPKHRYTP LKDNWVNIFT PIVKNLGLQI RFNLKKRQVE LRNPADREDT TDLQKATDFV RAFILGFEVN DAIALIRLDH LFLETFEIAD VKHSLKGDHV SRAIGRIAGK DGRTKLVIEN TTKTRIVVAN TKIHILGAYQ NLKLARNAVC SLILGSNPSK VYGNLRNMAS RGAERL
Source	E.coli
Target Names	pno-1
Protein Names	Recommended name: RNA-binding protein pno-1
Expression Region	1-276
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.