



Recombinant Bacillus cereus Ribosomal protein L11 methyltransferase (prmA)

Product Code	CSB-MP733943BAAD
Abbreviation	prmA
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q634M9
Product Type	Recombinant Protein
Immunogen Species	Bacillus cereus (strain ZK / E33L)
Purity	≥85% (SDS-PAGE)
Sequence	MKWSEISIHT TEEAVEAVSH ILHEAGASGV AIEDPAELTK EREQQYGEIY ALNPDEYP AE GVLKAYFPQ TDSLHETIAG VKSSIDVLPS YDIEIGTGNI TVNEVNEEDW ATAWKKYYHP VQISDTFTIV PTWEEYTPSS PEEKIIELDP GMAFGTGTHP TTTMCIRALE KTVQPGDTII DVG TGSGLVLS IAAAKLGASS VQAYDLDPVA VESAEMNVRL NKTDVVSVG QNSLLEGIEG PVDLIVANLL AEIILLFPED AARVVKSGGL FITSGIIAAK EKVISEALEK AGFTIEEVLR MEDWVAIIAR NA
Source	Mammalian cell
Target Names	prmA
Protein Names	Recommended name: Ribosomal protein L11 methyltransferase Short name=L11 Mtase EC= 2.1.1.-
Expression Region	1-312
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.