



Recombinant Mouse Myocyte-specific enhancer factor 2D (Mef2d)

Product Code	CSB-YP733979MO
Abbreviation	Mef2d
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q63943
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	≥85% (SDS-PAGE)
Sequence	MGRKKIQIQR ITDERNRQVT FTKRKFGLMK KAYELSVLCD CEIALIIFNH SNKLFQYAST DMDKVLLKYT EYNEPHESRT NADIIETLRK KGFNGCDSPE PDGEDSLEQS PLEDKYRRA SEELDGLFRR YGSSVPAPNF AMPVTVVSN QSSMQFSNPS SSLVTPSLVT SSLTDPRLLS PQQPALQRNS VSPGLPQRPA SAGAMLGGDL NSANGACPSV VGNGYVSARA SPGLLPVANG NSLNKVIKPAK SPPPPTHNTQ LGAPSRKPDV RVITSQGGKG LMHHLTEDHL DLNNAQRLGV SQSTHSLTTP VVSVATPSLL SQGLPFSSMP TAYNTDYQLP SAELSSLPAF SSPAGLALGN VTAWQQPQPP QQPQPPQPPQ SQPQPPQPPQ QQPPQQPHL VPVLSNLIP GSPLPHVGAA LVTTHPHIS IKSEPVSPSR ERSPAPPPA VFPAARPEPG EGLSSPAGGS YETGDRDDGR GDFGPTLGLL RPAPEPEAEG SAVKRMRLDT WTLK
Source	Yeast
Target Names	Mef2d
Protein Names	Recommended name: Myocyte-specific enhancer factor 2D
Expression Region	1-514
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.