



Recombinant Bacillus cereus Aspartate--ammonia ligase (asnA)

Product Code	CSB-MP733989BAAD
Abbreviation	asnA
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q63D03
Product Type	Recombinant Protein
Immunogen Species	Bacillus cereus (strain ZK / E33L)
Purity	≥85% (SDS-PAGE)
Sequence	MYQSLMTVRE TQIAIKEVKT FFEDQLAKRL ELFRVSAPLF VTKKSGGLNDH LNGVERPIEF DMLHSGEELE IVHSLAKWKR FALHEYGYEA GEGLYTNMNA IRRDEELDAT HSIYVDQWDW EKIVQKEWRT VDYLQKTVLT IYGIFKDLED HLFEKYPFLG KYLPEEIVFV TSQELEDKYP ELTPKDREHA IAKEHGAVFI IGIGDALRSG EKHDGRAADY DDWKLNGDIL FWHPVLQSSF ELSSMGIRVD SKSLDEQLTK TGEDFKREYD FHKGILEDLV PLTIGGGIGQ SRMCMYFLRK AHIGEVQSSV WPDDLREACK KENIHLF
Source	Mammalian cell
Target Names	asnA
Protein Names	Recommended name: Aspartate--ammonia ligase EC= 6.3.1.1 Alternative name(s): Asparagine synthetase A
Expression Region	1-327
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.