



Recombinant *Bacillus cereus* Ribose import ATP-binding protein RbsA (rbsA)

Product Code	CSB-MP734002BAAD
Abbreviation	rbsA
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q63FX9
Product Type	Recombinant Protein
Immunogen Species	<i>Bacillus cereus</i> (strain ZK / E33L)
Purity	≥85% (SDS-PAGE)
Sequence	MRIEMKNISK AFNGNPVLKN AQFMIETGEV HALMGENGAG KSTLMKILTG VYKRDGGTIT IDGQERTFKN AKEAEEYGIA FIHQELNILP NLTVAENMFL GKELMYGKTG ILRTRQMNAI AQQQLAELGL HVKGAMLAGE LSVGQQQIIE IAKALMTNAS VIIMDEPTAA LTDREIETLF TVINKLRKEG VSFVYISHRM EEIFSICDAI TILRDGEYVG KRSIPETSFD EVVSMMVGRS IGERYPERNS QISDVIFEMR NGTKKGKFEN VSFQVRKGEI LGVAGLMGAG RTDIMKSIFG YEPLDSGQIF INGQEVKIDS PIDAIRQRIA FITEDRKSEG LVLDFSIREN LALPNLESLS KGSVLSNELE QQFTADMMKL LNVKASSGEQ AVKSLSGGNQ QKVVIKWLG IHPQLLILDE PTRGVDVGAK KEIYSIMNKL TEEGDAVIMV SSELPEVLGM SDRVLVIHEG KVG GILEKDE ASQESIMALA TGGE
Source	Mammalian cell
Target Names	rbsA
Protein Names	Recommended name: Ribose import ATP-binding protein RbsA EC= 3.6.3.17
Expression Region	1-494
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.