



Recombinant Burkholderia pseudomallei Peptide methionine sulfoxide reductase MsrB (msrB)

Product Code	CSB-BP734058BXW
Abbreviation	msrB
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q63V23
Product Type	Recombinant Protein
Immunogen Species	Burkholderia pseudomallei (strain K96243)
Purity	≥85% (SDS-PAGE)
Sequence	MSGDRDDPRY PYPKDDAELR RRLTPMQYEV TQHAATEPPF TGEYTDTEDA GIYHCVVCGT ALFESGAKFH SGCGWPSYFK PIDGEVIDEK MDYTHGMTRV EVR CNQCGAH LGHVFEDGPR DKTGLRYCIN SAALNFEAKP ERK
Source	Baculovirus
Target Names	msrB
Protein Names	Recommended name: Peptide methionine sulfoxide reductase MsrB EC=1.8.4.12 Alternative name(s): Peptide-methionine (R)-S-oxide reductase
Expression Region	1-143
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.