



# Recombinant Burkholderia pseudomallei Glycerol kinase (glpK)

|                          |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                   |
|--------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Code</b>      | CSB-EP734066BXW-B                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                 |
| <b>Abbreviation</b>      | glpK                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                              |
| <b>Storage</b>           | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.                                                                                                                                                                                                                                                                                                                |
| <b>Uniprot No.</b>       | Q63X50                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                            |
| <b>Product Type</b>      | Recombinant Protein                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                               |
| <b>Immunogen Species</b> | Burkholderia pseudomallei (strain K96243)                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                         |
| <b>Purity</b>            | ≥85% (SDS-PAGE)                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                   |
| <b>Sequence</b>          | MQDQYILALD QGTTSSRAML FDRQGNIVSI AQKEFEQIYP QPGWVEHDPQ<br>EIWSTQAGVA AEAVTRVGLN GTAIAAIGIT NQRETTIVWD RETGHPINYA<br>IVWQDRRTAD FCDKLKAQGL SEKVRKTGL PIDSYFSATK IRWILDNVEG<br>AREKARQGKL AFGTVDSWLW WNFTKHELHV TDVTNASRTM LFNIHTLDWD<br>DELLDALEIP RSMLPQVRAS SEIYGPTKTT VFASKIPLAG IAGDQQAALF<br>GQMCTSSGMV KNTYGTGCFL MMNTGETPIE SRNNLVTIA WQVDGKVSYA<br>LEGSIFIAGA VVQWLRDGLG IKSASEIEA LAAGVPHTDG VYLVPFAFAGL<br>GAPHWNARAR GSLFGVTRGT TSAHLARAAL DSIAYQSLDV LKAMEADSGI<br>RIGELRVDGG ASANNLLMQF QADLLGVDAV RPRITETTAL GAAYLAGLAI<br>GYWQNVDELH SQWQLERRFA PSMQAEQVGS CLAGWQRAVR<br>AAKAWADDTQ |
| <b>Source</b>            | E.coli                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                            |
| <b>Target Names</b>      | glpK                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                              |
| <b>Protein Names</b>     | Recommended name: Glycerol kinase EC= 2.7.1.30 Alternative name(s):<br>ATP:glycerol 3-phosphotransferase Glycerokinase Short name= GK                                                                                                                                                                                                                                                                                                                                                                                                                                                             |
| <b>Expression Region</b> | 1-500                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                             |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                               |
| <b>Tag Info</b>          | Tag type will be determined during the manufacturing process.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                     |
| <b>Protein Length</b>    | full length protein                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                               |
| <b>Reconstitution</b>    | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.                                                                                                                                                                                       |
| <b>Shelf Life</b>        | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.                                                                                                                                                                                                                                                                                                                                                                                                                                                        |



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.