



Recombinant Mouse Ras-related protein Rab-34 (Rab34)

Product Code	CSB-BP734076MO
Abbreviation	Rab34
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q64008
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	≥85% (SDS-PAGE)
Sequence	MNILAPVRRD RVLAE LPQCL KKEAALHVRK DFHPRVTCAC QEHR TGTVGF KISKVIVVGD LSVGKTCLIN RFCKDTFDKN YKATIGVDFE MERFEVLGVP FSLQLWDTAG QERFKCIAST YYRGAQAI III VFNLNDVASL EHTKQWL TDA LKENDPSNVL LFLVGSKKDL STPAQYSLME KDALKVAQEI KAEYWAVSSL TGENVREFFF RVAALTFEAN VLADVEKSGA RHIADVVRIN SDDKNLYLTA SKKKATCCP
Source	Baculovirus
Target Names	Rab34
Protein Names	Recommended name: Ras-related protein Rab-34 Alternative name(s): Ras-related homolog Ras-related protein Rab-39 Ras-related protein Rah
Expression Region	1-259
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.