



Recombinant Rat Ubiquitin-like modifier-activating enzyme ATG7 (Atg7)

Product Code	CSB-EP734090RA-B
Abbreviation	Atg7
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q641Y5
Product Type	Recombinant Protein
Immunogen Species	Rattus norvegicus (Rat)
Purity	>85% (SDS-PAGE)
Sequence	MGDPGLSKLQ FAPFNSALDV GFWHELTQKK LNEYRLDEAP KDIKGYYYNG DSAGLPTRLT LEFSAFDMSA PTPARCCPAM GTLHNTNTLE AFKTADKLL LEQSANEIWE AIKSGAALEN PMLLNKFLLL TFADLKKYHF YYWFCCPALC LPESIPLIRG PVGLDQRLSP KQIQALEHAY DDLCRTEGVT ALPYFLFKYD DDTVLVSLK HYSDFFQGQR TKLTVGVYDP CNLTQHPGWP LRNFLVLA AH RWSGSFQSVE VLCFRDRTMQ GARDVTHSII FEVKLPEMAF SPDCPKAVGW EKNQKGGMGP RMVNLSGCMD PKRLAESSVD LNLKLMCWRL VPTLDLDKVV SVKCLLLGAG TLGCNVARTL MGWGVHRVTF VDNAKISYSN PVRQPLYEFE DCLGGGKPKA LAAERLQKI FPGVNASGFN MSIPMPGHPV NFSDVTMEQA RRDVEQLEEL IDSHDVIFLL MDTRESRWLP TVIAASKRKL VINAALGFD T FVVMRHGLKK PKQQGAGDLC PSHLVAPADL GSSLFANIPG YKLGCYFCND VVAPGDSTRD RTLDQQCTVS RPGLAVIAGA LVELMVSVL QHPEGGYAIA SSSDDRMNEP P TSLGLVPHQ IRGFLSRFDN VLPVSLAFDK CTACSSKVL D QYEQEGFTFL AKVFNSSHSF LEDLTGLTLL HQETQAAEIW DMSDEETV
Source	E.coli
Target Names	Atg7
Protein Names	Recommended name: Ubiquitin-like modifier-activating enzyme ATG7 Alternative name(s): ATG12-activating enzyme E1 ATG7 Autophagy-related protein 7 Short name= APG7-like
Expression Region	1-698
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Target Details	This gene was identified based on homology to Pichia pastoris GSA7 and Saccharomyces cerevisiae APG7. In the yeast, the protein appears to be



required for fusion of peroxisomal and vacuolar membranes. The protein shows homology to the ATP-binding and catalytic sites of the E1 ubiquitin activating enzymes.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.