



Recombinant *Bacillus licheniformis* Peptide chain release factor 1 (prfA)

Product Code	CSB-YP734230BQU
Abbreviation	prfA
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q65DV1
Product Type	Recombinant Protein
Immunogen Species	<i>Bacillus licheniformis</i> (strain ATCC 14580 / DSM 13 / JCM 2505 / NBRC 12200 / NCIMB 9375 / NRRL NRS-1264 / Gibson 46)
Purity	≥85% (SDS-PAGE)
Sequence	MLDRLKSID RYEKLNELLS DPEVVNDPKK LREYSKEQSD LQETVEVYRR YRSASEQLSD AKAMLEEKLD SDMREMVKEE ISELQEEIES LTDQLKVLLI PKDPNDDKNV IMEIRGAAGG EEAALFAGNL YRMYSRYAEL QGWKTEVMEA NMTGTGGYKE IIFMINGNGA YSRLKYENGA HRVQRPETE SGGRIHTSTA TVACLPEAEE VEVDIHEKDI RVDTFASSGP GGQSVNTTMS AVRLTHLPTG VVVSCQDEKS QIKNKEKAMK VLRARIYDKF QQEAQAEYDQ NRKSAVGTGD RSERIRTYNF PQNRVTDHRI GLTIQKLDQI LEGKLDEVID ALIVEDQASK LQQAES
Source	Yeast
Target Names	prfA
Protein Names	Recommended name: Peptide chain release factor 1 Short name= RF-1
Expression Region	1-356
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.