



Recombinant *Bacillus licheniformis* DNA-directed RNA polymerase subunit alpha (rpoA)

Product Code	CSB-MP734293BQU
Abbreviation	rpoA
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q65P79
Product Type	Recombinant Protein
Immunogen Species	<i>Bacillus licheniformis</i> (strain ATCC 14580 / DSM 13 / JCM 2505 / NBRC 12200 / NCIMB 9375 / NRRL NRS-1264 / Gibson 46)
Purity	≥85% (SDS-PAGE)
Sequence	MIEIEKPKIE TVEISDDAKY GKFVVEPLER GYGTTLGNSL RRILLSSLPG AAVTSIQIDG VLHEFSTIEG VVEDVTIIL NIKKLALKIY SEEKTLEID VQGEVVTA A DITHDSDVEI LNPDLHIATL GQNASFRVRL TAQRGRGYTP ADANKRDDQP IGVIPIDSIF TPVSRVSYQV ENTRVGQITN YDKLTLDVWT DGSTGPKEAI ALGSKILTEH LNIFVGLTDE AQHAEIMVEK EEDQKEKVL MTIEELDLSV RSYNCLKRAG INTVQELANK TEEDMMKVRN LGRKSLEEVK AKLEELGLGL RKDD
Source	Mammalian cell
Target Names	rpoA
Protein Names	Recommended name: DNA-directed RNA polymerase subunit alpha Short name= RNAP subunit alpha EC= 2.7.7.6 Alternative name(s): RNA polymerase subunit alpha Transcriptase subunit alpha
Expression Region	1-314
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.