



Recombinant tRNA (Ile)-lysidine synthase (tilS)

Product Code	CSB-MP734371BUJ
Abbreviation	tilS
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q65ZY6
Product Type	Recombinant Protein
Immunogen Species	Borrelia baviensis (strain ATCC BAA-2496 / DSM 23469 / PBI) (Borrelia baviensis)
Purity	≥85% (SDS-PAGE)
Sequence	MHFLDDNIQI KIDKFYKKNS LNKNRVIVAF SGGADSTLL LNLKYYLSNN IVAFYFAHFI RTDNEQNKEI EHVKGFCGLY NIALQVKKCD IDIKSESVRL GVSIEELARK CRYNALENAL KENDANYIAL AHNENDQIET IIMRFFQGSF LDGLSGIPGV NKNIIRPLLE VSRLEIENFL SLNNISFFVD STNAQDLYLR NRVRNLLPS IEKIFKGYQK CLKRISEFSK EFVDYFEKDE FFPVEKGYKY YSFDLKTFLN FPKYLVFRLI FKILNSEGIV AKVSYKALNE AFKVEINRKK NNVLLKTNDF FLEKRHNKIN LIFKRDEKDY KPFDFILEVG KWHLSLSLGI LLKCLECNTS SVSRLKCCSY EFRYKFFKDR FKAKKFFSKF IRCNPIYLMML LALDNRLIGI IDLNTLNLVW SEKSILKKIS ISLIGLLKE
Source	Mammalian cell
Target Names	tilS
Protein Names	Recommended name: tRNA(Ile)-lysidine synthase EC= 6.3.4.19 Alternative name(s): tRNA(Ile)-2-lysyl-cytidine synthase tRNA(Ile)-lysidine synthetase
Expression Region	1-440
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.