



Recombinant Triosephosphate isomerase (tpiA)

Product Code	CSB-EP734393BUI
Abbreviation	tpiA
Storage	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>
Uniprot No.	Q662V9
Product Type	Recombinant Protein
Immunogen Species	Borrelia baviensis (strain ATCC BAA-2496 / DSM 23469 / PBi) (Borrelia baviensis)
Purity	≥85% (SDS-PAGE)
Sequence	MRKTFLAGNW KMHYTSTEAS IVAKKIATEV KTLKDDFVIM ITPPFTALSK VSECIKGSNI LLGAQNMSYM ESGARTSEIS PSMLEFGVE YVILGHSECR LYLGETDEII NKKILTGLKH PFKYLILCVG ETLNERDSEK TLDVVLNQVK KGLDGVSESD IKRIILAYEP IWAIGTGKTA TKEEAEEVHK AIRLEIMKLY SKSVSDNVII QYGGSVNSSN VKELMNEPDI DGALIGGASL KAESFLSIIN NVL
Source	E.coli
Target Names	tpiA
Protein Names	Recommended name: Triosephosphate isomerase Short name= TIM EC= 5.3.1.1 Alternative name(s): Triose-phosphate isomerase
Expression Region	1-253
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>