



Recombinant Ribosomal protein L11 methyltransferase (prmA)

Product Code	CSB-EP734410YAH-B
Abbreviation	prmA
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q665E3
Product Type	Recombinant Protein
Immunogen Species	Yersinia pseudotuberculosis serotype I (strain IP32953)
Purity	≥85% (SDS-PAGE)
Sequence	MPWIQLKLNT TGNQAESLGD VLVESGAVSV TFQDTHDNPV FEPLPGETRL WGDTDVIGLY DAETDMADV VV AMLECHPQIG KGFHIEKIEQL EDKDWEREWM DNFHMPMRFGE RLWICPSWRD VPDPTAVNVM LDPGLAFGTG THPTTALCLQ WLDSL DLNGK TLIDFGCGSG ILAIAALKLG AARAIGIDID PQAIQASRDN AQRNGVSERL ELYLAKDQPA ELSADV VVAN ILAGPLRELA LLISVLPTTG GHLGLSGVLA TQAAGVAQAY EDKFILDPVA EKEEWCRITG IKK
Source	E.coli
Target Names	prmA
Protein Names	Recommended name: Ribosomal protein L11 methyltransferase Short name= L11 Mtase EC= 2.1.1.-
Expression Region	1-293
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.