



Recombinant 3,4-dihydroxy-2-butanone 4-phosphate synthase (ribB)

Product Code	CSB-MP734413YAH
Abbreviation	ribB
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q665V6
Product Type	Recombinant Protein
Immunogen Species	Yersinia pseudotuberculosis serotype I (strain IP32953)
Purity	≥85% (SDS-PAGE)
Sequence	MNQTLTSDFG TPVERVERAI DALRNGRGVM VLDDESRENE GDMVFAAEAM TLEQMALTIR HGSGIVCLCI TDERRQQLDL PMMVTHNSSQ FQTAFTVTIE AAEGVTTGVS AADRLTTIRK AIADNAKPAD LNRPGHVFPL RGQPGGVLSR RGHTEASIDL ATLAGYKPAG VLCELTNDG SMAHAPEVIA FAKLHDMPVV TIDDLAAYLQ SRAKKAS
Source	Mammalian cell
Target Names	ribB
Protein Names	Recommended name: 3,4-dihydroxy-2-butanone 4-phosphate synthase Short name= DHBP synthase EC= 4.1.99.12
Expression Region	1-217
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.