



Recombinant Arginine N-succinyltransferase (astA)

Product Code	CSB-EP734460YAH-B
Abbreviation	astA
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q66B20
Product Type	Recombinant Protein
Immunogen Species	Yersinia pseudotuberculosis serotype I (strain IP32953)
Purity	≥85% (SDS-PAGE)
Sequence	MMKVRPVERR DLADIFELAG KTGVGMTSLP QNEQHLAARI ERALNTWQGS LDPGEQGYLE VLEDSEQQKV VGVSAIEVAV GLNDPWYNFR VGTLVHASKA LNVYKSVPTL FLSNDHTGYS ELCTFLDPD YRKDKNGPFL SKVRFLFIAA FRQYFSRKVI AEMRGYTDEQ GRSPFWESVG RHFFSIEFAK ADYLSGTGQK AFIAELMPKH PLYVDFLAEE ARAVIGQVHP HTAPARAVLE TEGLYQGYV DIFDGGPTLE ANTDDVRAVR DSSKRTTVIK DYDIEDYDID PNGRLYLVA DHYHHFRAIL MNTHLSDERL RLTPESAEAL GVAAGDSVRI VSLFAPETKR
Source	E.coli
Target Names	astA
Protein Names	Recommended name: Arginine N-succinyltransferase Short name= AST EC= 2.3.1.109 Alternative name(s): AOST
Expression Region	1-350
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.