



Recombinant Probable endonuclease 4 (nfo)

Product Code	CSB-MP734470YAH
Abbreviation	nfo
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q66CR5
Product Type	Recombinant Protein
Immunogen Species	Yersinia pseudotuberculosis serotype I (strain IP32953)
Purity	≥85% (SDS-PAGE)
Sequence	MKFVGAHVSA AGGVDQAVIR AHELEATAFA LFTKNQRQWR AAPLAEDVIE KFKLACEKYG YTSAQILPHD SYLINLGHPV TEALEKSREA FIDELVRCQQ LGLSLLNFHP GSHLLQIDED QCLARIAESI NIALDATEGV TAVIENTAGQ GSNLGFKFEH LAIIERVED KSRVGVICIDT CHAFAAGYDL RTEEDCEHTF AALGKIVGFQ YLRGMHLNDA KSEFNSRVDR HHSLGEGNIG KTVFSYIMRD SRFDNIPLIL ETVNMDIWAE EIAWLKSQAE IEPSL
Source	Mammalian cell
Target Names	nfo
Protein Names	Recommended name: Probable endonuclease 4 EC= 3.1.21.2 Alternative name(s): Endodeoxyribonuclease IV Endonuclease IV
Expression Region	1-285
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.