



Recombinant Probable endonuclease 4 (nfo)

| | |
|--------------------------|---|
| Product Code | CSB-YP734470YAH |
| Abbreviation | nfo |
| Storage | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |
| Uniprot No. | Q66CR5 |
| Product Type | Recombinant Protein |
| Immunogen Species | Yersinia pseudotuberculosis serotype I (strain IP32953) |
| Purity | ≥85% (SDS-PAGE) |
| Sequence | MKFVGAHVSA AGGVDQAVIR AHELEATAFA LFTKNQRQWR AAPLAEDVIE KFKLACEKYG YTSAQILPHD SYLINLGHPV TEALEKSREA FIDELVRCQQ LGLSLLNFHP GSHLLQIDED QCLARIAESI NIALDATEGV TAVIENTAGQ GSNLGFKFEH LAIIERVED KSRVGVICIDT CHAFAAGYDL RTEEDCEHTF AALGKIVGFQ YLRGMHLNDA KSEFNSRVDR HHSLGEGNIG KTVFSYIMRD SRFDNIPLIL ETVNMDIWAE EIAWLKSQAE IEPSL |
| Source | Yeast |
| Target Names | nfo |
| Protein Names | Recommended name: Probable endonuclease 4 EC= 3.1.21.2 Alternative name(s): Endodeoxyribonuclease IV Endonuclease IV |
| Expression Region | 1-285 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | full length protein |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| Shelf Life | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |