



Recombinant Rat NEDD4-binding protein 2-like 2 (N4bp2l2)

Product Code	CSB-BP734515RA
Abbreviation	N4bp2l2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q66H65
Product Type	Recombinant Protein
Immunogen Species	Rattus norvegicus (Rat)
Purity	≥85% (SDS-PAGE)
Sequence	MPYSEIEAKF LGPGKELTRE PCYKKLKS AV DDGVFPHRGG PDSHRIQ EKT KNNRLPVATI NFRRRVCPGE DKTSTDVLKP LHKEMPGDKV GGTESIGSQA LQDGKPLAPA RDDEIYSTSK AFIGPIYKPP EKKKCRERKN ETATFNNTDS KRRQEEKQKF NSKKLEIDTE LSQFYKEIEE LENENEASQG SCKEPEPSQE QTLSHDQACN TLKPEGESKD LSTVLQSHCG YQPYLEDES D YPCDEQLIPA FCDTSFSPFR PEWQSVHPFV IPHDPLPSFN YFNIQRFGAP LHPPPSVFHA RDDARLQSGC YIDSYQDGWS CLTLDQNDEY ANYDVPGNVH PFGNGYGVQD DSVKNGFCEI RECWQDPSMD MHNEADRFVN QCFQEDRLNK LQKLLILLRG LPGSGKTTLS RILLGQSRDG IVFSTDDYFH HQDGYRYNVN QLGDAHDWNQ NRAKQAIDQG RSPVIIDNTN TQAWEMKPYV EMAIGKGYRV EFHEPETWWK FDPEELEKRN KHGVSRRKIA QMLDRYEFQM SISIVMNSVE PTQKSTQTPL PLQGDQRWGG SLGSHSQVSI TDDY
Source	Baculovirus
Target Names	N4bp2l2
Protein Names	Recommended name: NEDD4-binding protein 2-like 2 Alternative name(s): Phosphonoformate immuno-associated protein 5 homolog
Expression Region	1-574
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.