



# Recombinant Rat NEDD4-binding protein 2-like 2 (N4bp2l2)

|                          |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                               |
|--------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Code</b>      | CSB-MP734515RA                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                |
| <b>Abbreviation</b>      | N4bp2l2                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                       |
| <b>Storage</b>           | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.                                                                                                                                                                                                                                                                                                                                                                                                            |
| <b>Uniprot No.</b>       | Q66H65                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                        |
| <b>Product Type</b>      | Recombinant Protein                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                           |
| <b>Immunogen Species</b> | Rattus norvegicus (Rat)                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                       |
| <b>Purity</b>            | ≥85% (SDS-PAGE)                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                               |
| <b>Sequence</b>          | MPYSEIEAKF LGPGKELTRE PCYKKLKS AV DDGVFPHRGG PDSHRIQ EKT<br>KNNRLPVATI NFRRRVCPGE DKTSTDVLKP LHKEMPGDKV GGTESIGSQA<br>LQDGKPLAPA RDDEIYSTSK AFIGPIYKPP EKKKCRERKN ETATFNNTDS<br>KRRQEEKQKF NSKKLEIDTE LSQFYKEIEE LENENEASQG SCKEPEPSQE<br>QTLSHDQACN TLKPEGESKD LSTVLQSHCG YQPYLEDESD YPCDEQLIPA<br>FCDTSFSPFR PEWQSVHPFV IPHDPLPSFN YFNIQRFGAP LHPPPSVFHA<br>RDDARLQSGC YIDSYQDGWS CLTLDQNDEY ANYDVPGNVH<br>PFGNGYGVQD DSVKNGFCEI RECWQDPSMD MHNEADRFVN<br>QCFQEDRLNK LQKLLILLRG LPGSGKTTLS RILLGQSRDG IVFSTDDYFH<br>HQDGYRYNVN QLGDAHDWNQ NRAKQAIDQG RSPVIIDNTN<br>TQAWEMKPYV EMAIGKGYRV EFHEPETWWK FDPEELEKRN KHGVSRRKIA<br>QMLDRYEFQM SISIVMNSVE PTQKSTQTPL PLQGDQRWGG SLGSHSQVSI<br>TDDY |
| <b>Source</b>            | Mammalian cell                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                |
| <b>Target Names</b>      | N4bp2l2                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                       |
| <b>Protein Names</b>     | Recommended name: NEDD4-binding protein 2-like 2 Alternative name(s):<br>Phosphonoformate immuno-associated protein 5 homolog                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                 |
| <b>Expression Region</b> | 1-574                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                         |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                           |
| <b>Tag Info</b>          | Tag type will be determined during the manufacturing process.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                 |
| <b>Protein Length</b>    | full length protein                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                           |
| <b>Reconstitution</b>    | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.                                                                                                                                                                                                                                                                                   |



## Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.