



Recombinant Rat Protein lyl-1 (Lyl1)

Product Code	CSB-EP734524RA
Abbreviation	Lyl1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q66HH3
Product Type	Recombinant Protein
Immunogen Species	Rattus norvegicus (Rat)
Purity	≥85% (SDS-PAGE)
Sequence	MCPPQAQAEV GSAMTEKTEM VCTSSPAPAP PSKPASPGPL PAEEVDHRNT CNPWLPPGVP VINLGHTRPT GAAMPTTELS AFRPSLLQLA ALGTAPPTLA LHYHHPFLN SVYIGPAGPF SIFPNSRLKR RPSHGELDLV DGHQPQKVAR RVFTNSRERW RQQHVNGAFA ELRKLPTHP PDRKLSKNEV LRLAMKYIGF LVRLLRDQAA VLASGPSAPG SRKPPAHRGV EGNARCGAGH RVEAARSQPV LPGDCDGDPN GSVRPIKMEQ TALSPEVR
Source	E.coli
Target Names	Lyl1
Protein Names	Recommended name: Protein lyl-1 Alternative name(s): Lymphoblastic leukemia-derived sequence 1
Expression Region	1-278
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.