



Recombinant Danio rerio Poly (ADP-ribose) glycohydrolase ARH3 (adprhl2)

Product Code	CSB-MP734525DIL
Abbreviation	adprhl2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q66HT8
Product Type	Recombinant Protein
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Purity	≥85% (SDS-PAGE)
Sequence	MSAAARVVAP VMMLSRFRGA LVGSVLGDCI GGEFEGAVDV PLDRVLQHLS ALEDDTRGDG ILQYSDDTAM MRCVADSLLT RMTFDERDMA QRFAKEYSHS PGRGYGSGVV QVLRKLASPH LKDVFPQAQA QFGGRGSFGN GGAMRAVPFA LAFRSRADVR KYSRFGAMLT HSCSLGYNGA ALQALAVHLS LQGALALPKD FIDKLISEME ELEKDETAKH DAKALNLSEF PYCSRLHRVK ELMDKTSVSI EEVISELGNG IAALQSVPTA IFCVLYCLEP QDGLPERFGG LERTIAYSLA LGGDTDTIAC MAGAIAGAHY GIDSIPLSWQ VSCEGVDEAD DLARRLYDLY CLPQHNE DRG NNQPHTTNTD
Source	Mammalian cell
Target Names	adprhl2
Protein Names	Recommended name: Poly(ADP-ribose) glycohydrolase ARH3 EC= 3.2.1.143 Alternative name(s): ADP-ribosylhydrolase 3 [Protein ADP-ribosylarginine] hydrolase-like protein 2
Expression Region	1-370
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.