



Recombinant Trachypithecus obscurus Mannose-binding protein C (MBL2)

Product Code	CSB-EP734549TKL
Abbreviation	MBL2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q66S41
Product Type	Recombinant Protein
Immunogen Species	Trachypithecus obscurus (Dusky leaf-monkey) (Presbytis obscura)
Purity	>85% (SDS-PAGE)
Sequence	ETVTCEDSQK ICPAVIACNS PGINGFPGKD GRDGTKGEKG EPGQGLRGLQ GPPGKLGPPG NPGPSGSPGP KGQKGDGPES PDCDSSLAAS ERKALQTEMA HIKKWLTFSL GRQVGNKFFL TNGEMMTFDK VKALCAKFQA SVATPRNAE NRAIQNLIKE EAFLGITDEN TEGQFVDLIG NRLTYTNWNN GEPNAGSDE DCVLLLKNGK WNDVPCSSSH LALCEFPI
Source	E.coli
Target Names	MBL2
Protein Names	Recommended name: Mannose-binding protein C Short name= MBP-C Alternative name(s): MBP1 Mannan-binding protein Mannose-binding lectin
Expression Region	21-248
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.