



Recombinant Influenza A virus Matrix protein 1 (M)

Product Code	CSB-BP734557IMJ
Abbreviation	M
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q67173
Product Type	Recombinant Protein
Immunogen Species	Influenza A virus (strain A/Equine/Tennessee/5/1986 H3N8)
Purity	≥85% (SDS-PAGE)
Sequence	MSLLTEVETY VLSIVPSGPL KAEIAQRLED VFAGKNTDLE ALMEWLKTRP ILSPLTKGIL GFVFTLTVPS ERGLQRRRFV QNALSGNGDP NNMDRAVKLY RKLKREITFH GAKEVALSYS TGALASCMGL IYNRMGTVTT EVAFGLVCAT CEQIADSQHR SHRQMVTTTN PLIRHENRMV LASTTAKAME QMAGSSEQAA EAMEVASKAR QMVQAMRTIG THPSSSAGLK DDLENLQAY QKRMGVQMQR FK
Source	Baculovirus
Target Names	M
Protein Names	Recommended name: Matrix protein 1 Short name= M1
Expression Region	1-252
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.