



# Recombinant Adenylosuccinate synthetase (purA)

|                          |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                    |
|--------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Code</b>      | CSB-BP734587SRW                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                    |
| <b>Abbreviation</b>      | purA                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                               |
| <b>Storage</b>           | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.                                                                                                                                                                                                                                 |
| <b>Uniprot No.</b>       | Q67J55                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                             |
| <b>Product Type</b>      | Recombinant Protein                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                |
| <b>Immunogen Species</b> | Symbiobacterium thermophilum (strain T / IAM 14863)                                                                                                                                                                                                                                                                                                                                                                                                                                                                |
| <b>Purity</b>            | ≥85% (SDS-PAGE)                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                    |
| <b>Sequence</b>          | MPVVVVMGAQ WGDEGKGKQV DLLAERAQVV VRSTGGSNAG<br>HTVWAGGRQY KLHQVPSGIL YPGTLCVIGH GVVLDPPKLL EEMDRLASQG<br>VDLSSLRISG GAHIVFPFHI RLDEAEEDRK GDRKIGTTRR GIGPAYMDKF<br>ARVGIRLVDM LDRDEFLPKL TALVEEKNRI LEKVYGLPGF TVEEIAEPYL<br>EYAEERLRPYV ANTVELVNDA IDAGKNVLF GAQGHLLDID FGTPYVTAS<br>HPIAAGAIIG AGVGPTKVS R VGVVKAYTS RVGEGPFPTE LHGEEAHRIR<br>EEGHEYGTTT GRPRRIGWLD LVMVRYACRV SGITDLAVPH LDTLAKTGLP<br>TLKVCVGYRM PDGTVTREFP VGLKALSQVE PVYEELPNWE WSAEMSTARQ<br>YDELPEGARR YVRLIEDVTG VRVSILGVGS ERTQAIYRSA IF |
| <b>Source</b>            | Baculovirus                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                        |
| <b>Target Names</b>      | purA                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                               |
| <b>Protein Names</b>     | Recommended name: Adenylosuccinate synthetase Short name= AMPSase<br>Short name= AdSS EC= 6.3.4.4 Alternative name(s): IMP--aspartate ligase                                                                                                                                                                                                                                                                                                                                                                       |
| <b>Expression Region</b> | 1-432                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                              |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.                                                                                                                                                                                                                                                                                                                                                                                                                |
| <b>Tag Info</b>          | Tag type will be determined during the manufacturing process.                                                                                                                                                                                                                                                                                                                                                                                                                                                      |
| <b>Protein Length</b>    | full length protein                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                |
| <b>Reconstitution</b>    | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.                                                                                                        |
| <b>Shelf Life</b>        | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.                                                                                                                                                                                                                                 |