



Recombinant Adenylosuccinate synthetase (purA)

Product Code	CSB-EP734587SRW
Abbreviation	purA
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q67J55
Product Type	Recombinant Protein
Immunogen Species	Symbiobacterium thermophilum (strain T / IAM 14863)
Purity	≥85% (SDS-PAGE)
Sequence	MPVVVVMGAQ WGDEGKGKQV DLLAERAQVV VRSTGGSNAG HTVWAGGRQY KLHQVPSGIL YPGTLCVIGH GVVLDPPKLL EEMDRLASQG VDLSSLRISG GAHIVFPFHI RLDEAEEDRK GDRKIGTTRR GIGPAYMDKF ARVGIRLVDM LDRDEFLPKL TALVEEKNRI LEKVYGLPGF TVEEIAEPYL EYAEERLPYV ANTVELVNDI IDAGKNVLF GAQGHLLDID FGTYPYVTAS HPIAAGAIIG AGVGPTKVSF VVGVVKAYTS RVGEGPFPTE LHGEEAHRIR EEGHEYGTTT GRPRRIGWLD LVMVRYACRV SGITDLAVPH LDTLAKTGLP TLKVCVGYRM PDGTVTREFP VGLKALSQVE PVYEELPNWE WSAEMSTARQ YDELPEGARR YVRLIEDVTG VRVSILGVGS ERTQAIYRSA IF
Source	E.coli
Target Names	purA
Protein Names	Recommended name: Adenylosuccinate synthetase Short name= AMPSase Short name= AdSS EC= 6.3.4.4 Alternative name(s): IMP--aspartate ligase
Expression Region	1-432
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.