



Recombinant Alanine racemase (alr)

Product Code	CSB-BP734596SRW
Abbreviation	alr
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q67K79
Product Type	Recombinant Protein
Immunogen Species	Symbiobacterium thermophilum (strain T / IAM 14863)
Purity	≥85% (SDS-PAGE)
Sequence	MDLDGMRPTW AEVDLNAVRA NIRALKRISK APRLMAVVKA NGYGHGAVPV ATAAIEAGAD WLGVASVEEG VTLRRHGISA PILVLGYVSP GQAEAVLTEG LRVALFDGEL GQALNREGRR LGRWARVHLK VDTGMGRIGL QPAEAVGRLGR ELARLDHVEV EGVFTHLATA DEPGNPYTRL QLERYEAALA ELAAAGVRPA IRHAANSAGL MLHPEAHYDM VRSGIAVVGL PPAPGVAWPV KLAPALTWKT RVGLVKWLEA GHSISYGCTY TTARREQUIAT LPVGYADGYP RRLSNRAQVL IRGRRCPVVG VVTMDQMMVR VPDDL PVRVG DEVVLIGRQG GEEITATELA GLADTISYEI VCGISRRVPR FYGGETA
Source	Baculovirus
Target Names	alr
Protein Names	Recommended name: Alanine racemase EC= 5.1.1.1
Expression Region	1-377
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.