



Recombinant Chaperone protein DnaJ (dnaJ)

Product Code	CSB-BP734624SRW
Abbreviation	dnaJ
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q67S53
Product Type	Recombinant Protein
Immunogen Species	Symbiobacterium thermophilum (strain T / IAM 14863)
Purity	≥85% (SDS-PAGE)
Sequence	MDRLAKRDYY EILGVPRNAT EAEIKKAFRN LARKYHPDAN KDDPDAAEFK KEINEAYQVL SDPEKRARYD QFGHAAEQMG GAGGNPFEGF GGFGDFGGFS DIFEMFFGGA GRQRNPRGPV RGADLEYELE LTLKEAAFGC TKELRVPRVE DCDTCHGSGA RPGTQPVTCP KCGGTGQVQM TQHTVGRFV NVMTCDRCRG EGKIVESPCP TCRGRGRVQK TQRVEVKVPG GVETGTRLRM PGYGEAGERG GPPGDLFIVM RVRPDRRFRR EGDDLFTTTE ISFIQAALGT EIEVETLDGP ELIKIPEGTQ PGDTIRLKGK GTHRLRGSGR GDLHVVISVK TPGRLSERER ELLLEVAALR GERVAGMTEN KEKSFLKKMK DALGGR
Source	Baculovirus
Target Names	dnaJ
Protein Names	Recommended name: Chaperone protein DnaJ
Expression Region	1-386
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.