



Recombinant ATP synthase subunit beta (atpD)

Product Code	CSB-BP734630SRW
Abbreviation	atpD
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q67TB7
Product Type	Recombinant Protein
Immunogen Species	Symbiobacterium thermophilum (strain T / IAM 14863)
Purity	≥85% (SDS-PAGE)
Sequence	MSEPTGKILA IMGAVVDVEF PDGHLPKLNH ALVIRNGGER GDINLTLEVA QHLGNNQVRC IAMDSTDGLR RGETVYDTQA PITVPVGEPV LGRIFNVLGE TIDGKGPVNS PVSLPIHREP PGVANVDPAT EMLPTGIKVI DLLCPYARGG KIGLFGGAGV GKTVMIQELI HNIAYKFGGY SIFAGVGERT REGNDLYHEM IEAGVIDKVA MVFGQMNEPP GARMRVALTG LTMAEYFRDE QGLDVLLFID NIFRFTQAGS EVSALLGRMP SAVGYQPTLA TEMGYLQERI TSTKKGSITS VQAIYVPADD LTD PAPATTF AHL DATTVLS RAVFEKAIFP AVDPLDSTSR ILDPNIVGRE HYEVARGVQK VLQRYKELQD IAILGMDEL SDEDKLIVAR ARKIERFLSQ PMFVAEKFTN LPGVFVDPKD TVKGFKEILE GKHDDLPEQA FYMVGTIEQA VEKGGKLLAE VS
Source	Baculovirus
Target Names	atpD
Protein Names	Recommended name: ATP synthase subunit beta EC= 3.6.3.14 Alternative name(s): ATP synthase F1 sector subunit beta F-ATPase subunit beta
Expression Region	1-472
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.