



Recombinant Arabidopsis thaliana

Uncharacterized protein At4g14100 (At4g14100)

Product Code	CSB-YP734638DOA
Abbreviation	At4g14100
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q67YC9
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	DEPVP TPWPHQFHAL MFMNYSGLDLS MIDLWYDWTN GRNFNIIQEQ LGGITYDLEW NNGTSFFYTL DESKSCRSQG LEVGILRPNW LDGAKYLGQQ NVSGFLCNVW EKVDFIWYYE DVETKRPVQW IFYTGREAHI MTYEVGAVLE DEKWQAPVYC FNKEKSLST KGALSEFKGY KEKAAL
Source	Yeast
Target Names	At4g14100
Protein Names	Recommended name: Uncharacterized protein At4g14100
Expression Region	26-206
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.