



Recombinant Human Olfactomedin-like protein 2B (OLFML2B)

Product Code	CSB-YP734650HU
Abbreviation	OLFML2B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q68BL8
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	IVLTGTSE PPDAQTVAPA EDETLQNEAD NQENVLSQLL GDYDKVKAMS EGSDCQCKCV VRPLGRDACQ RINAGASRKE DFYTVETITS GSSCKCACVA PPSALNPCEG DFRLQKLREA DSQDLKLSTI IDMLEGAFYG LDLLKLHSVT TKLVGRVDKL EEEVSKNLTK ENEQIKEDME EIRTEMNKRG KENCSENILD SMPDIRSALQ RDAAYAHYHP EYEERFLQEE TVSQQINSIE LLQTRPLALP EYVKSQRPLQ RQVHLRGRPA SQPTVIRGIT YYKAKVSEEE NDIEEQQDEF FSGDNGVDLL IEDQLLRHNG LMTSVTRRPA ATRQGHSTAV TSDLNARTAP WSSALPQPST SDPSIANHAS VGPTLQTTSV SPDPTRESVL QPSPQVPATT VAHTATQQPA APAPPAVSPR EALMEAMHTV PVPPTTVRTD SLGKDAPAGW GTTTPASPTLS PEEEDDIRNV IGRCKDTLST ITGPTTQNTY GRNEGAWMKD PLAKDERIYV TNYYYGNTLV EFRNLENFKQ GRWSNSYKLP YSWIGTGHVV YNGAFYYNRA FTRNIIKYDL KQRYVAAWAM LHDVAYEEAT PWRWQGHSDV DFAVDENGLW LIYPALDDEG FSQEVIVLSK LNAADLSTQK ETTWRTGLRR NFYGNCFVIC GVLYAVDSYN QRNANISYAF DTHNTQIVP RLLFENEYSY TTQIDYNPKD RLLYAWDNGH QVTYHVIFAY
Source	Yeast
Target Names	OLFML2B
Protein Names	Recommended name: Olfactomedin-like protein 2B Alternative name(s): Photomedin-2
Expression Region	23-750
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final



concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.