



# Recombinant Human HAUS augmin-like complex subunit 3 (HAUS3)

<b>Product Code</b>	CSB-BP734652HU
<b>Abbreviation</b>	HAUS3
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q68CZ6
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	SCGNEFVET LKKIGYPKAD NLNGEDFDWL FEGVEDESFL KWFCGNVNEQ NVLSERELEA FSILQKSGKP ILEGAALDEA LKTCKTSDLK TPRLDDKELE KLEDEVQTLK KLKLNKIQRN NKCQLMASVT SHKSLRLNAK EEEATKCLKQ SQGILNAMIT KISNELQALT DEVTQLMMFF RHSNLGQGTM PLVFLSQFSL EKYLSQEEQS TAALTLYTKK QFFQGIHEVV ESSNEDNFQL LDIQTPSICD NQEILEERRL EMARLQLAYI CAQHQLIHLK ASNSSMKSSI KWAEESLHSL TSKAVDKENL DAKISLTSE IMKLEKEVTQ IKDRSLPAVV RENAQLLNMP VVKGDFDLQI AKQDYITARQ ELVLNQLIKQ KASFELLQLS YEIELRKHRD IYRQLENLVQ ELSQSNMMLY KQLEMLTDPS VSQQINPRNT IDTKDYSTHR LYQVLEGENK KKELFLTHGN LEEVAEKLKQ NISLVQDQLA VSAQEHSFFL SKRNKDVDMML CDTLYQGGNQ LLLSDQELTE QFHKVESQLN KLNHLLTDIL ADVKTKRRTL ANNKLHQMER EFVYFLKDE DYLKDIVENL ETQSKIKAVS LED
<b>Source</b>	Baculovirus
<b>Target Names</b>	HAUS3
<b>Protein Names</b>	Recommended name: HAUS augmin-like complex subunit 3
<b>Expression Region</b>	2-603
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



## Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.