



Recombinant Human Coiled-coil domain-containing protein 102B (CCDC102B)

Product Code	CSB-EP734653HU-B
Abbreviation	CCDC102B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q68D86
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MNLDSIHRLI EETQIFQMQQ SSIKSRGDMV APASPPRDTC NTCFPLHGLQ SHAAHNFCAH SYNTNKWDIC EELRLRELEE VKARAAQMEK TMRWWSICTA NWREKWSKVR AERNSAREEG RQLRIKLEMA MKELSTLKKK QSLPPQKEAL EAKVTQDLKL PGFVEESCEH TDQFQLSSQM HESIREYLVK RQFSTKEDTN NKEQGVVIDS LKLSEEMKPN LDGVDLFNNG GSGNGETKTG LRLKAINLPL ENEVTEISAL QVHLDEFQKI LWKEREMRTA LEKEIERLES ALSLWKWKYE ELKESKPKNV KEFDILLGQH NDEMQELSGN IKEESKSQNS KDRVICELRA ELERLQAENT SEWDKREILE REKQGLEREN RRLKIQVKEM EELLDKKNRL SANSQSPDFK MSQIDLQEKN QELLNLQHAY YKLNRYQAN IAELTHANNR VQNEAEVKK LRLRVEELKQ GLNQKEDELD DSLNRQIRKLQ RSLDEEKERN ENLETELRLH QNW
Source	E.coli
Target Names	CCDC102B
Protein Names	Recommended name: Coiled-coil domain-containing protein 102B
Expression Region	1-513
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.